CE
0483

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Luxor Z True Nature

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(EN) Instruction for use

Luxor Z True Nature

Please read these instructions for use carefully before using the product.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use. The product must only be used by dentists, dental technicians and adequately trained dental staff. Original tools and components must be used for processing.

1. Product description / Intended purpose

Luxor Z True Nature zirconium oxide blanks for the manufacture of dentures are manufactured using CIP technology and presintered at low temperatures. The blanks are available in Vita classical A1 - D4 and bleach. Although the product has a certain strength, please handle it with care. When the product is received, check it for completeness or damage. In the event of any defects, please contact the responsible customer advisor or the bredent contact address provided.

2. Indication

Individual abutments, veneers, inlays and onlays, anterior and posterior crowns and bridges, occlusally screwed crowns and bridges, monolithic crowns and bridges, circular bridges.

3. Contraindication

Based on the current state of knowledge there are no known contraindications. The product must not be used in the event of a known or suspected incompatibility of one of its components.

4. Hazard and precautionary statements

The batch numbers of all products used must be recorded to ensure traceability and for the purpose of handling complaints.

Hazard statements:

Not currently known

4.1 Prevention

Read and understand all the safety information and product-accompanying documents before use.

Use personal protective equipment as required. Protective equipment includes:

- Protective clothing
- Protective gloves
- Mouth-nose protection, respiratory protection
- Protective goggles and shields

5. Storage and shelf life information

Storage:

Dry storage.

For the purpose of identification and protection of the product, the product must be stored in its original outer packaging (delivery condition) (no storage of units that have already been unpacked or partially unpacked).

Shelf life:

The product has no shelf life limitation.

6TH Application

6.1 Preparation guidelines (guide)

Luxor Z True Nature is not suitable for off-the-shoulder or tangential preparations as both methods result in a thinner margin and risk of fracture.

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The inner edge of shoulder must be slick or be fluted.

Preparation width of incisal, labial and palatal should be above 0,7 mm.

The axial wall must be blunt (milling radius) in order to be able to achieve a suitable milling result with CAD/CAM.

Monolithic:

The inner edge of shoulder must be slick or be fluted.

The minimum width of incisal edge, labial side, and palatal side is at least 0,7 mm.

Occlusal surface width should be circular min. 1,0 mm.

The axial wall must be blunt in order to achieve an ideal milling effect using CAD/CAM.

6.2 Scanning

Oral Scanning:

Before oral scanning, the first thing is to do gingival treatment to exposure the preparation area and clean and dry the blood and saliva in scanning area.

First scan the preparation area, then scan the adjacent teeth or the entire jaw.

Note: when scanning the edge, the speed should be slow and repeat scanning.

Scanning the abutment:

Scan each tooth position step by step, ensuring that the teeth are in full contact with the model.

6.3 Design

For monolithic restorations, certain requirements must be met as follows:

Restoration type	Anterior tooth		Posterior tooth	
	Thickness (minimum)	Connectors	Thickness (minimum)	Connectors
Single crown	0,7 mm	-	0,7 mm	-
3-unit bridge	0,8 mm	9 mm ²	1,0 mm	12 mm ²
Bridges with 4 or more units	1,2 mm	12 mm ²	1,2 mm	12 mm ²

1. The zirconium bridge design ensures crack resistance, so the height of the connectors should be as high as possible. The cross-section of the connector for a 3-unit anterior bridge should be at least 9 mm² and 12 mm² for a 4-unit anterior bridge. The cross-section of the posterior connector should be at least 12 mm² .

2. A maximum of 2 pontics between the abutment teeth/implants.

3. Use cutter radial correction. The upper part of the abutment should be filled with wax before performing the scan and design.

6.4 Nesting

Select zirconium oxide blank with a suitable thickness.

The connectors should be placed on the buccal side in the most protruding position. The thickness of the connectors should be at least 2 mm. In order to achieve a better effect of tooth colour and translucency, it is recommended to have 0,5 mm space between the top and bottom layer of the restoration and the blank, e.g. the restoration et the ébauche, pour une restauration qui par ex. serait de 13,7 mm – then a blank with a height of 16 mm should be chosen. In the case of bridges, it must be ensured that the bridge connectors are not positioned in the upper third.


6.5 Inserting the blank

Please note the blank labelling, which contains important information:

- Material type
- Height
- Colour

Please note the label on the side of the blank, which contains important information:

- LOT
- Directional arrow pointing in the incisal direction (the side with the highest translucency)
- Scaling factor

6.6 Milling

When restoring with Luxor Z True Nature, make sure that an undamaged bre.cam Cutter ZR or DCS red.line tool is used and the zirconium oxide is only processed dry. After completing the milling, check for the defects listed below:

- Are there any cracks?
- Are there any impurities?
- Is anything broken?

If one of these defects occurs, the reasons for them must be identified and the restorations milled again.

6.7 Removing the restoration

Use a technician's handpiece with suitable grinding tools to separate restorations from the blanks. A towel should be placed on the table before grinding to prevent restorations from falling, snapping, or breaking on the table.

Perforate all connectors 90% without pressure, in order to avoid tilting the restorations, then completely cut through the connectors without pressure. However, it is not recommended to make too many adjustments to the soft restorations to avoid hidden cracks or chipping etc.

6.8 Cleaning

Thoroughly remove zirconium dust on the surface and inside of restorations with a ceramic brush.

CAUTION
The following effects may occur if cleaning/advice is not observed: 1. After sintering, the milling residue/zirconium dust adheres to the surface and inside the restorations. In addition to the white dots, which may change the aesthetics, the fit of the restorations is impaired.

6.9 Sintering

Step 1: Checking the sintered beads

CAUTION
1. The beads must be replaced in the following situations: strong colour changes, incomplete or damaged shape, dull colour, lack of light transmission.
2. Check that the beads stick together, if so, separate the sintered beads to ensure good mobility.
3. Quantity: Completely cover the bottom of the tray with 2-3 layers.
4. If new beads are used, sinter the beads with residual material 1-2 times and then sinter the restorations.

Step 2: Checking the sintering furnace

Cleaning and maintenance of the furnace

1. Cleaning method: Remove impurities in the combustion chamber and place in the furnace with discarded zirconium oxide residues for sintering (cleaning fire). Clean once a week.

2. Check heating elements for defects before sintering. If there are any chippings on the surface of the element (MoSi2), place the separated zirconium oxide residues in the furnace for cleaning firing.

3. When the furnace is not in use, keep it closed to ensure a dry furnace interior.

4. Keep sintering equipment safe and clean, avoid dust formation as this can affect the heating elements.

Strictly follow the sintering parameters as listed below.

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