Ivotion Dent/Ivotion Dent Multi

[en]	Instructions for Use DCL tooth material for the CAD/CAM technique	[hr]	Upute za uporabu DCL materijal zuba za CAD/CAM tehniku
[de]	Gebrauchsinformation DCL-Zahnmaterial für die CAD/CAM-Technik	[cs]	Návod k použití Dentální DCL materiál pro techniku CAD/CAM
[fr]	Mode d'emploi Matériau dentaire DCL pour la technique CAD/CAM	[sk]	Návod na používanie Zubný materiál DCL pre techniku CAD/CAM
[it]	Istruzioni d'uso Materiale dentale DCL per la tecnica CAD/CAM	[hu]	Használati utasítás DCL foganyag CAD/CAM technikához
[es]	Instrucciones de uso Material para dientes DCL mediante la técnica CAD/CAM	[sr]	Упутство за употребу DCL материјал за зубе за CAD/CAM технику
[pt]	Instruções de Uso Material odontológico DCL para técnica CAD/CAM	[mk]	Упатство за употреба DCL забен материјал за CAD/CAM техника
[sv]	Bruksanvisning DCL tandmaterial för CAD/CAM-tekniken	[bg]	Инструкции за употреба DCL материал за зъби за CAD/CAM техника
[da]	Brugsanvisning DCL-tandmateriale til CAD/CAM-teknikken	[sq]	Udhëzimet e përdorimit Material dhëmbi DCL për teknikën CAD/CAM
[fi]	Käyttöohjeet DCL-materiaali CAD/CAM-menetelmään	[ro]	Instrucțiuni de utilizare Material dentar DCL pentru tehnica CAD/CAM
[no]	Bruksanvisning DCL-tannmateriale til CAD/CAM-teknologi	[uk]	Інструкція з використання Зубний матеріал DCL (поперечно зшитий полімер) для систем CAD/CAM
[nl]	Gebruiksaanwijzing DCL-tandmateriaal voor de CAD/CAM-techniek	[et]	Kasutusjuhend DCL-i hambamaterjal CAD-/CAM-tehnika jaoks
[el]	Οδηγίες Χρήσης Υλικό τεχνητών δοντιών DCL για την τεχνική CAD/CAM	[lv]	Lietošanas instrukcija DCL zobu materiāls CAD/CAM metodei
[tr]	Kullanım Talimatları CAD/CAM tekniği için DCL diş materyali	[lt]	Naudojimo instrukcija DCL dantų medžiaga, skirta CAD / CAM technikai
[ru]	Инструкция по применению Материал для зубов DCL (поперечно-сшитый полимер) для использования по технологии CAD/CAM		
[pl]	Instrukcja stosowania Materiał protetyczny DCL do techniki CAD/CAM		
[sl]	Navodila za uporabo Zobni material DCL za tehniko CAD/CAM		

CE 0123

Rx ONLY

Date information prepared:
2024-01-23 / Rev. 0

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

ivoclar

1 Intended use

Intended purpose

Fabrication of artificial teeth for removable dentures

Patient target group

- Patients with permanent teeth
- Adult patients with dental implants
- Edentulous adult patients

Intended users

- Dental laboratory technicians (fabrication of restorations in the dental laboratory)
- Denturists (fabrication of restorations in the dental laboratory)

Special training

No special training required.

Use

For dental use only.

Description

Ivotion Dent and Ivotion Dent Multi are tooth-coloured discs for the CAD/CAM technique to fabricate denture teeth for removable denture prosthetics using a subtractive manufacturing process. Ivotion Dent consists of one layer, and Ivotion Dent Multi consists of three layers. The discs are made of highly cross-linked polymethyl methacrylate (PMMA), which is also called DCL (Double Cross-Linked) material. DCL material consists of a modified PMMA variant, in which both the previously produced polymer and the matrix are highly cross-linked. The discs are available in a variety of shades to faithfully reproduce the individual characteristics of the patient.

Technical data

	According to ISO 10477:2018*	Typical materials properties
Flexural strength	≥ 80 MPa	120 MPa
Water absorption	≤ 40 µg/mm ³	25 µg/mm ³
Water solubility	≤ 7.5 µg/mm ³	0.8 µg/mm ³

* Currently, there is no standard that defines the requirements for machinable tooth materials processed with milling systems. The product complies with the requirements of ISO 22112:2017. The physical data are determined according to ISO 10477:2018.

Indications

- Partial edentulism in the anterior and posterior region
- Complete edentulism

Types of restorations:

- Complete dentures
- Partial dentures
- Implant-supported dentures
- Tooth-supported dentures

Contraindications

The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Observe the required minimum thickness (2.5 mm).
- Only mill the discs in a compatible CAD/CAM system.
- Staining / layering with SR Nexco (light-curing, MMA-based conditioner for the bonding of light-curing materials to PMMA materials (heat-curing or self-curing polymers and denture teeth)).
- Avoid overheating of the material when finishing the objects.
- Due to its high translucency, the material is not suitable for the fabrication of crowns, bridges and hybrid abutments for esthetic reasons.
- Avoid contact with solvents.

- Protect the discs from direct sunlight once you have removed them from the packaging.
- Failure to observe the stipulated limitations of use may lead to failure.
- The final denture is for use by one patient only.

Side effects

In individual cases, allergic reactions to methyl methacrylate materials have been reported.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Highly cross-linked PMMA, pigments

2 Application

NOTICE! Avoid overheating the material through improper use of the milling tools, polishing brushes and steam cleaners.

NOTICE! Ivotion Dent and Ivotion Dent Multi are based on the A-D shade system.

Designing the denture

NOTICE! Observe the following specifications:

- Wall thickness for teeth: min 2.5 mm

NOTICE! The position of the teeth within the Ivotion Dent Multi disc determines the esthetics of the restoration.

1. Design the denture.

Fabricating the teeth

NOTICE! Use a compatible milling system with a standard holder of 98.5 mm.

NOTICE! Only use milling tools recommended for the particular milling system.

NOTICE! Make sure that the disc holder is securely fixed.

NOTICE! Ivotion Dent Multi must be placed in the milling machine with the printed side pointing towards the incisal/occlusal region of the restoration.

1. Mill the teeth.

Preparing the teeth

1. Detach the holding bars. Use a suitable tungsten carbide bur.
2. Smooth out the attachment points.

Optional: Finishing the teeth

NOTICE! When finishing segments manually, make sure that the connector dimension remains as large as possible.

1. Finish the teeth. Use a cross-cut or diamond tungsten carbide bur.

Checking the fit

1. Check the fit of the teeth in the denture base.
2. In case of interfering contacts: Adjust the tooth. Use a cross-cut or diamond tungsten carbide burs.

Bonding the teeth to the denture base

1. Bond the teeth on the denture base. Use a suitable bonding system (e.g. Ivotion Bond).  For guidance regarding processing, consult the Instruction for Use of the product.

Finishing/polishing the denture

1. Finish the denture. Use cross-cut or diamond tungsten-carbide burs.

NOTICE! Gloss varnish must not be used.

2. Optional: Characterize the denture. Use the lab composite SR Nexco.  For guidance regarding processing, consult the Instruction for Use of the product.

NOTICE! When polishing the teeth, preserve the anatomical surface structure.

3. Polish the denture. Use a goat-hair brush and pumice.
4. Polish the denture to a high gloss. Use a cotton or leather polishing buff and polishing paste.

Cleaning the denture

1. Clean the denture under running water. Use a tooth brush.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use and the explanation of symbols are available on the website: www.ivoclar.com/elFU
- Basic UDI-DI: 76152082APROS008W9

Warnings

- Observe the Safety Data Sheet (SDS) (available at Ivoclar Vivadent AG www.ivoclar.com).
- Do not inhale grinding dust. Use an extraction unit and wear a face mask.
- Wear protective goggles during finishing.
- Processing the material may result in sharp edges. There is a risk of injury.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

None

4 Shelf life and storage

- Protect the product from direct sunlight.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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- Perskaitykite saugos duomenų lapą (SDL) (jį rasite Ivoclar Vivadent AG www.ivoclar.com).